




NASOLABIAL INJECTION



MECHANISM OF ACTION



Singfiller® shows equivalent results on Suitability, Safety & Effectivity.



Data from: "Summary of Clinical Investigation for Cross-linked Sodium Hyaluronate Gel for Plastic Surgery"



Singfiller®

Cross-linked sodium Hyaluronate for Plastic Surgery



www.singclean.net
www.singfiller.org

Hangzhou Singclean Medical Products Co., Ltd
 Add: No. 125 (E), 10th Street, Hangzhou Economic and
 Technological Development Zone, Zhejiang, China
 P.C.: 310018
 Tel: +86-571-63431868 Fax: +86-571-63431886
 E-mail: sales@hzxhe.com
 Web: www.singclean.net www.singfiller.org

Singfiller®

sterile injection
 safety and natural
 strong resistance


CE
 1023



COMPANY
PROFILE

Hangzhou Singclean Medical Products Co., Ltd established in March 2003 with registered capital of USD13,000,000 is a national high-tech enterprise integrating the R&D, production and sales of medical biomaterials. With well-qualified staff, 55 Standard and first-class facilities, the company ensures the continuous production of medical devices products of Class III that confirms to the laws and regulations of the State Food and Drug Administration (SFDA) and EU MDD.

Sticking to the quality policy of “Quality First, Management Priority, Building a Harmonious Enterprise, Creating the Singclean Brand”, the company passed the ISO9001 and ISO13485 quality management system in 2007, and obtained the Innovation Fund Certificate by Ministry of Science and Technology with the Microporous Polysaccharide Hemospheres project in 2008. The company obtained the certification of China high-tech enterprise in 2010. The medical sodium hyaluronate gels passed the CE certificate in 2011, and the cross-linked sodium hyaluronate gel for plastic surgery brand— **Singfiller** ® passed the CE certificate in 2017.



WHAT IS SINGFILLER?

Singfiller is BDDE-cross-linked hyaluronic acid filler which is colorless, transparent, viscous, sterile non-pyrogenic injection. This third generation product is reticulated and obtained via bio-fermentation technology.

Singfiller is developed with cutting-edge technology, guaranteed for purity, safety and long-lasting effect.



TYPE	Fine	Derm	Deep	Deeper	Sub-Q
COMPOSITION	20mg/ml(HA)	20mg/ml(HA)	20mg/ml(HA)	20mg/ml(HA)	20mg/ml(HA)
Particle size	0.10~0.15mm	0.15~0.28mm	0.28~0.50mm	0.50~1.00mm	1.25~2.00mm
Recommended indications	thin superficial lines ⑦⑧⑨⑩	moderate wrinkle ①④⑤⑥⑦⑧	Deep facial wrinkles and folds ①⑥⑦⑧⑨⑩	Deep facial wrinkles and folds ①⑥⑧	the local depression,soft tissue filler. And breast enhance ,breast butt deep dermis or subcutaneous tissue shallow
Where to inject	upper part of dermis	middle part of dermis	deep layer of dermis	deep dermis or subcutaneous tissue shallow	
Volume of Syringe	1mL (CE)	1mL (CE) / 2mL (CE)	1mL (CE) / 2mL (CE)	1mL (CE) / 2mL (CE) / 10mL(No CE)	2mL(No CE) / 20mL(No CE)
Storage	0~30℃ (32-86°F)	0~30℃ (32-86°F)	0~30℃ (32-86°F)	0~30℃ (32-86°F)	0~30℃ (32-86°F)
Use	Single	Single	Single	Single	Single
Duration	6 - 9months	6 - 9months	6 - 9months	6-12months	6-12months
Needle	2.25PRTC	27G*1	26G*1	26G*1	/
	1.0PRTC	30G*2	27G*2	26G*2	/

APPLICATION

⚠️ facial beauty injection

⚠️ improve skin appearance

⚠️ structure,reduce wrinkle,prevent skin aging

⚠️ improve labial shape

⚠️ facial tissue augmentation

